

# Healthy Products From Nose To Tail!™

Provide Your Four-Legged Customers Quality and Value at a Price Your Pet Parents Can Afford!



## Made Specially for Dogs and Cats!™

The Healthy Tracks for Pets™ Portfolio and Test Buddy™ Pet-Monitoring Blood Glucose System are specially designed to fit the needs of dogs and cats with diabetes. At an affordable price, our high-quality<sup>1</sup> products offer accurate results and reliable performance to help dogs and cats with diabetes live a happier & healthier life.

Order today! For more product information, contact your Trividia Health, Inc. representative at 1-800-877-6161.

Healthy Tracks™  
for Pets

TEST BUDDY™

TRIVIDIA  
HEALTH™



## Test Buddy™ Pet-Monitoring Blood Glucose System

The Test Buddy™ Pet-Monitoring Blood Glucose System was specially designed with a specific chemistry and unique algorithm for pet parents to track and monitor their pet's blood glucose levels confidently. Test Buddy™ APP allows pet parents to sync glucose results from the Test Buddy™ meter memory to certain mobile devices.



They can easily track, monitor and share their pet's blood glucose results with their veterinarian.

The Test Buddy™ Pet-Monitoring Blood Glucose Test Strips are only for use with the Test Buddy™ Pet-Monitoring Meter, and feature a beveled tip, and



two fill-detect electrodes to ensure adequate blood sample and prevent false starts and erroneous results.

### Features

- No coding
- About 10 seconds
- Stores 1000 results with date/time
- Small 0.5 microliter sample size
- 7-, 14-, 30-, 60- and 90-day averaging

- Ear or Paw Pad Testing Sites



- Event tagging



## Healthy Tracks for Pets™ Lancing Device and Lancets

The Healthy Tracks for Pets™ Lancing Device, provides an easy and efficient way for pet parents to draw a successful blood sample from their dog or cat with diabetes. The lancing device is for use with the Healthy Tracks for Pets™ Lancets, available in different sizes and designed with a tri-beveled needle tip to make the experience more comfortable.



### Lancing Device Features

- Easy to use
- 5 adjustable depth settings
- Ergonomic, tri-lobular design for comfortable grip

### Lancets Features

- Simple-to-use twist top for easy testing
- Sterile, single use
- Available in 26G and 28G



## Healthy Tracks for Pets™ Syringes

The Healthy Tracks for Pets™ Syringes are designed to administer insulin to dogs and cats with diabetes. These single-use syringes come in a variety of gauges and sizes to help pet parents deliver their pet's insulin with ease.

### Features

- Large clear unit markings for accurate dosing
- Sterile, single use
- Non-toxic, non-pyrogenic and latex free
- Available in 29G and 31G
- U-40 and U-100 syringes available

Test Buddy™ Pet-Monitoring Blood Glucose System	UPC
Meter Kit	0-21292-01365-3
Test Strips, 50 ct.	0-21292-01369-1
Health Tracks for Pets™	
Lancing Device	0-21292-01380-6
Sterile Lancets, 100ct., 26G	0-21292-01367-7
Sterile Lancets, 100ct., 28G	0-21292-01379-0
Syringes U-40 29G 1/2" 0.3cc	0-21292-01368-4
Syringes U-40 29G 1/2" 0.5cc	0-21292-01399-8
Syringes U-40 29G 1/2" 1.0cc	0-21292-01400-1
Syringes U-40 29G 1/2" 0.3cc half unit	0-21292-01401-8
Syringes U-40 29G 1/2" 0.5cc half unit	0-21292-01402-5
Syringes U-100 31G 0.3cc	0-21292-01406-3
Syringes U-100 31G 0.5cc	0-21292-01407-0
Syringes U-100 31G 1.0cc	0-21292-01408-7

Visit [www.healthy-tracks.com](http://www.healthy-tracks.com) for more information.